

Boule Con-Diff

CONTROL

IVD



Control	Lot	
Tri-Level	22401-1K	2024-05-30
Low	22401-11	2024-05-30
Normal	22401-12	2024-05-30
High	22401-13	2024-06-04

Open vial stability 14 days ¹

MEDONIC		CONTROL	L
Medonic M-series/M32		LOT	22401-11
		Mean ²	± Range ³
RBC	10 ¹² /L	2.23	0.12
MCV	fL	66.1	5.0
RDW	%	14.2	4.0
Hct	%	14.7	2.5
Plt	10 ⁹ /L	79	15
MPV	fL	11.1	1.5
WBC	10 ⁹ /L	3.3	0.3
Hgb	g/dL	5.6	0.3
MCH	pg	25.1	1.5
MCHC	g/dL	38.1	3.0
Lymph#	10 ⁹ /L	1.5	0.4
Mid#	10 ⁹ /L	0.2	0.3
Gran#	10 ⁹ /L	1.6	0.8
Lymph	%	45.8	7.0
Mid	%	5.0	5.0
Gran	%	49.2	7.0

CONTROL	N
LOT	22401-12
Mean ²	± Range ³
4.18	0.18
78.5	5.0
14.2	4.0
32.8	3.0
225	30
11.4	1.5
8.3	0.6
11.7	0.4
28.0	2.5
35.7	3.0
3.9	0.6
0.4	0.5
4.0	1.0
46.7	5.0
5.0	5.0
48.3	5.0

CONTROL	H
LOT	22401-13
Mean ²	± Range ³
5.08	0.22
91.7	5.0
13.6	4.0
46.6	3.5
513	60
11.9	1.5
20.6	1.8
16.2	0.5
31.9	1.5
34.8	3.0
12.7	1.5
1.0	1.0
6.9	2.5
61.5	5.0
5.0	5.0
33.5	5.0

- 1
- bg Трайност на отворен флакон 14 дни (bg)
 - hr Stabilnost otvorene bočice 14 dani (hr)
 - cs Stabilita otevřené ampulky 14 dní (cs)
 - da Åben rør stabilitet 14 dage (da)
 - et Stabiilsus avatud vialis 14 päevade (et)
 - fr Stabilité en flacon ouvert 14 jour (fr)
 - de Stabilität geöffneter Flaschen 14 tage (de)
 - el Διατήρηση δείγματος μετά την αποσφράγιση 14 ημέρα (el)
 - it Stabilità della fiala aperta 14 giorni (it)
 - lv Atvērt flakonu stabilitāti 14 dienas (lv)
 - lt Stabilumas atidarius buteliuką 14 dienas (lt)
 - no Åpen rør/glass stabilitet 14 dager (no)
 - pl Trwałość po otwarciu fiołki 14 Liczba dni (pl)
 - pt Estabilidade após abertura do frasco 14 dias (pt)
 - ro Termenul de valabilitate al fiolei desfăcute 14 zile (ro)
 - ru Стабильность вскрытого флакона 14 дней (ru)
 - sr Proizvod stabilan nakon otvaranja pakovanja 14 dana (sr)
 - sk Stabilita otvorenej injekčnej liekovky 14 dní (sk)
 - sl Stabilita otvorenej injekčnej liekovky 14 dni (sl)
 - es Estabilidad de la cápsula abierta 14 días (es)
 - sv Hållbarhet för öppnad flaska 14 dagar (sv)
 - tr Açık şişe stabilitesi 14 gün (tr)

- 2
- Средно
 - Sredstvo
 - Průměr
 - Middelverdi
 - Keskmine
 - Moyenne
 - Mittelwert
 - Μέση
 - Media
 - Nozīmē
 - Vidurkis
 - Middelverdi
 - Średnia
 - Média
 - Medie
 - Среднее
 - Srednja vrednost
 - Priemer
 - Povprečno
 - Media
 - Medelvärde
 - Ortalama

- 3
- Обхват
 - Raspon
 - Rozsah
 - Område
 - Vahemik
 - Intervalle
 - Bereich
 - Αναμενόμενο εύρος α
 - Intervallo
 - Diapazons
 - Intervalas
 - Område
 - Zakres
 - Intervalo
 - Interval
 - диапазон
 - Opseg
 - Rozsah
 - Območje
 - Intervalo
 - Intervall
 - Aralığı

For Medonic M-series / M32 with barcode reader

M-Series:

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

M32:

1. In the Main Menu go to "Quality Control" "Input Assays"
2. Scan in each blood calibrator / control column by column, sequence 1 to 9

Instr.: M-series
Lot: 22401-11
Prod.: Con-Diff Low
ID: 2240111+
Exp.: 30-May-2024

Instr.: M-series
Lot: 22401-12
Prod.: Con-Diff Norm
ID: 2240112+
Exp.: 30-May-2024

Instr.: M-series
Lot: 22401-13
Prod.: Con-Diff High
ID: 2240113+
Exp.: 4-Jun-2024

