

# Boule Con-Diff

CONTROL

IVD



Open vial stability 14 days <sup>1</sup>

Control	Lot	
Tri-Level	22404-1K	2024-08-29
Low	22404-11	2024-08-29
Normal	22404-12	2024-08-29
High	22404-13	2024-09-02

MEDONIC		CONTROL	L
Medonic M-series/M32		LOT	22404-11
		Mean <sup>2</sup>	± Range <sup>3</sup>
RBC	10 <sup>12</sup> /L	2.21	0.12
MCV	fL	67.5	5.0
RDW	%	14.4	4.0
Hct	%	14.9	2.5
Plt	10 <sup>9</sup> /L	84	15
MPV	fL	11.6	1.5
WBC	10 <sup>9</sup> /L	3.3	0.3
Hgb	g/dL	5.6	0.3
MCH	pg	25.3	1.5
MCHC	g/dL	37.6	3.0
Lymph#	10 <sup>9</sup> /L	1.5	0.4
Mid#	10 <sup>9</sup> /L	0.2	0.3
Gran#	10 <sup>9</sup> /L	1.6	0.8
Lymph	%	45.0	7.0
Mid	%	5.8	5.0
Gran	%	49.2	7.0

CONTROL	N
LOT	22404-12
Mean <sup>2</sup>	± Range <sup>3</sup>
4.19	0.18
81.3	5.0
14.8	4.0
34.1	3.0
229	30
12.0	1.5
8.2	0.6
11.9	0.4
28.4	1.5
34.9	3.0
3.7	0.6
0.6	0.5
3.9	1.0
44.8	5.0
7.6	5.0
47.6	5.0

CONTROL	H
LOT	22404-13
Mean <sup>2</sup>	± Range <sup>3</sup>
5.11	0.22
92.9	5.0
14.1	4.0
47.5	3.5
517	60
12.0	1.5
20.5	1.8
16.2	0.5
31.7	1.5
34.1	3.0
12.4	1.5
1.2	1.0
6.9	2.5
60.5	5.0
5.7	5.0
33.8	5.0

- 1
- bg Трайност на отворен флакон 14 дни (bg)
  - hr Stabilnost otvorene bočice 14 dni (hr)
  - cs Stabilita otevřené ampulky 14 dní (cs)
  - da Åben rør stabilitet 14 dage (da)
  - et Stabiilsus avatud viaalis 14 päevade (et)
  - fr Stabilité en flacon ouvert 14 jour (fr)
  - de Stabilität geöffneter Flaschen 14 tage (de)
  - el Διατήρηση δείγματος μετά την αεροσφράγιση 14 ημέρα (el)
  - it Stabilità della fiala aperta 14 giorni (it)
  - lv Atvērt flakonū stabilitāti 14 dienas (lv)
  - lt Stabilumas atidarius buteliuką 14 dienas (lt)
  - no Åpen rør/glass stabilitet 14 dager (no)
  - pl Trwałość po otwarciu fiołki 14 Liczba dni (pl)
  - pt Estabilidade após abertura do frasco 14 dias (pt)
  - ro Termenul de valabilitate al fiolei desfăcute 14 zile (ro)
  - ru Стабильность вскрытого флакона 14 дней (ru)
  - sr Proizvod stabilan nakon otvaranja pakovanja 14 dana (sr)
  - sk Stabilita otvorenej injekčnej liekovky 14 dní (sk)
  - sl Stabilita otvorenej injekčnej liekovky 14 dni (sl)
  - es Estabilidad de la cápsula abierta 14 días (es)
  - sv Hållbarhet för öppnad flaska 14 dagar (sv)
  - tr Açık şişe stabilitesi 14 gün (tr)

- 2
- Sredno
  - Sredstvo
  - Průměr
  - Middelværdi
  - Keskmine
  - Moyenne
  - Mittelwert
  - Μέση
  - Media
  - Nozīmē
  - Vidurkis
  - Middelverdi
  - Šrednia
  - Média
  - Medie
  - Среднее
  - Srednja vrednost
  - Priemer
  - Povprečno
  - Media
  - Medelvärde
  - Ortalama

- 3
- Обхват
  - Raspon
  - Rozsah
  - Område
  - Vahemik
  - Intervalle
  - Bereich
  - Αναμενόμενο εύρος α
  - Intervallo
  - Diapazons
  - Intervalas
  - Område
  - Zakres
  - Intervalo
  - Interval
  - диапазон
  - Opseg
  - Rozsah
  - Območje
  - Intervalo
  - Intervall
  - Aralığı

# For Medonic M-series / M32 with barcode reader

## M-Series:

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

## M32:

1. In the Main Menu go to "Quality Control" "Input Assays"
2. Scan in each blood calibrator / control column by column, sequence 1 to 9

Instr.:	M-series	Instr.:	M-series	Instr.:	M-series
Lot:	22404-11	Lot:	22404-12	Lot:	22404-13
Prod.:	Con-Diff Low	Prod.:	Con-Diff Norm	Prod.:	Con-Diff High
ID:	2240411+	ID:	2240412+	ID:	2240413+
Exp.:	29-Aug-2024	Exp.:	29-Aug-2024	Exp.:	2-Sep-2024

1:	 Bou l e C B A V 2	 Bou l e C B A V 2	 Bou l e C B A V 2
2:	 2240411+	 2240412+	 2240413+
3:	 24082900001153030000	 24082900002153030000	 24090200003153030000
4:	 10004030800000000000	 10004030000000000000	 10004030800000000000
5:	 22112675500560314925	 41918813501190434130	 51122929501620547535
6:	 08415116152531537630	 22930120152841534930	 51760120153171534130
7:	 03303015040020301608	 08206037060060503910	 20518124150121006925
8:	 14440450700585049270	 14840448500765047650	 14140605500575033850
9:	 11117317337558331174	 11111312030812249474	 11112465960227665154