

Boule Con-Diff

CONTROL

IVD



Control	Lot	
Tri-Level	22407-1K	2024-11-28
Low	22407-11	2024-11-28
Normal	22407-12	2024-11-28
High	22407-13	2024-12-02

Open vial stability 14 days ¹

SWELAB		CONTROL	L
Swelab Alfa/Alfa Plus		LOT	22407-11
		Mean ²	± Range ³
RBC	10 ¹² /L	2.26	0.12
MCV	fL	67.6	5.0
RDW	%	14.4	4.0
Hct	%	15.3	2.5
Plt	10 ⁹ /L	85	15
MPV	fL	12.1	1.5
WBC	10 ⁹ /L	3.2	0.3
Hgb	g/dL	5.6	0.3
MCH	pg	24.8	1.5
MCHC	g/dL	36.6	3.0
Lymph#	10 ⁹ /L	1.6	0.4
Mid#	10 ⁹ /L	0.1	0.3
Gran#	10 ⁹ /L	1.5	0.8
Lymph	%	49.2	7.0
Mid	%	4.1	5.0
Gran	%	46.7	7.0

CONTROL	N	
LOT	22407-12	
	Mean ²	± Range ³
	4.23	0.18
	86.2	5.0
	13.2	4.0
	36.5	3.0
	245	30
	12.5	1.5
	8.4	0.6
	12.6	0.4
	29.8	1.5
	34.5	3.0
	4.3	0.6
	0.4	0.5
	3.7	1.0
	51.2	5.0
	4.4	5.0
	44.4	5.0

CONTROL	H	
LOT	22407-13	
	Mean ²	± Range ³
	5.14	0.22
	92.9	5.0
	12.8	4.0
	47.8	3.5
	525	60
	12.7	1.5
	21.1	1.8
	16.1	0.5
	31.3	1.5
	33.7	3.0
	13.5	1.5
	0.8	1.0
	6.8	2.5
	63.9	5.0
	3.9	5.0
	32.2	5.0

- 1
- bg Трайност на отворен флакон 14 дни (bg)
 - hr Stabilnost otvorene bočice 14 dani (hr)
 - cs Stabilita otevřené ampulky 14 dní (cs)
 - da Åben rør stabilitet 14 dage (da)
 - et Stabiilsus avatud viaalis 14 päevade (et)
 - fr Stabilité en flacon ouvert 14 jour (fr)
 - de Stabilität geöffneter Flaschen 14 tage (de)
 - el Διατήρηση δείγματος μετά την αποσφράγιση 14 ημέρα (el)
 - it Stabilità della fiala aperta 14 giorni (it)
 - lv Atvērt flakonu stabilitāti 14 dienas (lv)
 - lt Stabilumas atidarius buteliuką 14 dienas (lt)
 - no Åpen rør/glass stabilitet 14 dager (no)
 - pl Trwałość po otwarciu fiołki 14 Liczba dni (pl)
 - pt Estabilidade após abertura do frasco 14 dias (pt)
 - ro Termenul de valabilitate al fiolei desfăcute 14 zile (ro)
 - ru Стабильность вскрытого флакона 14 дней (ru)
 - sr Proizvod stabilan nakon otvaranja pakovanja 14 dana (sr)
 - sk Stabilita otvorenej injekčnej liekovky 14 dní (sk)
 - sl Stabilita otvorenej injekčneje liekovky 14 dni (sl)
 - es Estabilidad de la cápsula abierta 14 días (es)
 - sv Hållbarhet för öppnad flaska 14 dagar (sv)
 - tr Açık şişe stabilitesi 14 gün (tr)

- 2
- Средно
 - Sredstvo
 - Průměr
 - Middelverdidi
 - Keskmine
 - Moyenne
 - Mittelwert
 - Μέση
 - Media
 - Nozīmē
 - Vidurkis
 - Middelverdi
 - Średnia
 - Média
 - Medie
 - Среднее
 - Srednja vrednost
 - Priemer
 - Povprečno
 - Media
 - Medelvärde
 - Ortalama

- 3
- Обхват
 - Raspon
 - Rozsah
 - Område
 - Vahemik
 - Intervalle
 - Bereich
 - Αναμενόμενο εύρος α
 - Intervallo
 - Diapazons
 - Intervalas
 - Område
 - Zakres
 - Intervalo
 - Interval
 - диапазон
 - Orseg
 - Rozsah
 - Območje
 - Intervalo
 - Intervall
 - Aralığı

For Swelab Alfa / Alfa Plus with barcode reader

Alfa:

1. In the Main Menu go to "Q/C Menu" "Enter Con/Cal" "Control/Calibrator Assay Value Input"
2. Scan in each blood calibrator/control column by column, sequence 1 to 9

Alfa Plus:

1. In the Main Menu go to "Quality Control" "Input Assays"
2. Scan in each blood calibrator / control column by column, sequence 1 to 9

Instr.: Swelab Alfa
Lot: 22407-11
Prod.: Con-Diff Low
ID: 2240711+
Exp.: 28-Nov-2024

Instr.: Swelab Alfa
Lot: 22407-12
Prod.: Con-Diff Norm
ID: 2240712+
Exp.: 28-Nov-2024

Instr.: Swelab Alfa
Lot: 22407-13
Prod.: Con-Diff High
ID: 2240713+
Exp.: 2-Dec-2024

